



Product Datasheet

IRAK4 antibody (orb689092)

Description	Rabbit monoclonal antibody to IRAK4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	A synthesized peptide derived from human IRAK4
Target	IRAK4
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q9NWZ3
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Jurkat cells...



Western blot
analysis of
THP-1 whole
cel...