


## Product Datasheet

### S6 Ribosomal Protein (phospho-Ser240/244) antibody (orb6891)

## Description

S6 Ribosomal Protein (phospho-



## Species/Host

Rabbit

## Reactivity

Human, Mouse, Rat

## Conjugation

Unconjugated

## Tested

## Applications

FC, IHC-P, WB

Immunohistochemical  
staining of human  
lu...

## Immunogen

KLH conjugated Synthesised  
phosphopeptide derived from human RPS6  
around the phosphorylation site of  
Ser240/244 (RA(p-S)TSK(p-S)ES)

## Target

Phospho-RPS6 (Ser240 + Ser244)

## Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin,  
0.03% Proclin300 and 50% Glycerol.

## Form/Appearance

Liquid

## Concentration

1mg/ml

## Storage

Shipped at 4°C. Store at -20 °C for one  
year. Avoid repeated freeze/thaw cycles.

## Note

For research use only

## Isotype

IgG

## Clonality

Polyclonal

## MW

29 kDa

## Uniprot ID


**P62753**

## Dilution Range

WB=1:100-1000 , IHC-P=1:100-500, IHC-  
F=1:100-500, IF=1:100-500, ELISA=1:100-  
1000

## Expiration Date

12 months from date of receipt.



Immunohistochemical  
staining of rat brai...