


Product Datasheet

ADAM17 antibody (orb689125)

Description	Rabbit monoclonal antibody to ADAM17
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	A synthesized peptide derived from human ADAM17
Target	ADAM17
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P78536
Dilution Range	Recommended dilution: IF:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of HeLa
Cell...