

Product Datasheet

SARS Antibody (orb689356)

Description	Anti-SARS-CoV-2 Nucleocapsid antibody(DM37), Rabbit mAb
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Nucleocapsid
Form/Appearance	Liquid
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	Rabbit IgG
Clonality	Recombinant
Uniprot ID	PODTC9
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.