



Product Datasheet

RSK2 (phospho-Tyr529) antibody (orb6897)


Description	RSK2 (phospho-Tyr529) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RSK2 around the phosphorylation site of Tyr529 (VE(p-Y)LH)
Target	Phospho-RSK2 (Tyr529)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	84 kDa
Uniprot ID	P51812
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of A673 Cell US...



Flow cytometric analysis of (mo) Brain C...



Immunohistochemical staining of human ki...