



Product Datasheet

S100A6 Antibody (orb68992)

Description	Mouse monoclonal antibody to S100A6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Purified recombinant fragment of calyculin expressed in E. Coli.
Preservatives	Purified antibody in PBS containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Entrez	6277
Dilution Range	IHC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of truncated calyculin...



Immunohistochemical analysis of paraffin-embedded tissue...