

Product Datasheet

RSK3 (phospho-Thr356/Ser360) antibody (orb6900)

Description	RSK3 (phospho-Thr356/Ser360) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RSK3 around the phosphorylation site of Thr356/Ser360 (AR(p-T)PTD(p-S)PG)
Target	Phospho-RSK3 (Thr356 + Ser360)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	83 kDa
Uniprot ID	Q15349
Expiration Date	12 months from date of receipt.