



Product Datasheet

SNCG (breast cancer-specific protein 1) Antibody (orb69000)

Description	Mouse monoclonal antibody to SNCG (breast cancer-specific protein 1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Purified recombinant fragment of SNCG expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Entrez	6623
Dilution Range	IHC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of truncated SNCG ...



Immunohistochemical analysis of paraffin...