



## Product Datasheet

### Myoglobin Antibody (orb69004)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Myoglobin                           |
| <b>Species/Host</b>        | Mouse                                                            |
| <b>Reactivity</b>          | Human                                                            |
| <b>Conjugation</b>         | Unconjugated                                                     |
| <b>Tested Applications</b> | ELISA, IHC                                                       |
| <b>Immunogen</b>           | Purified recombinant fragment of Myoglobin expressed in E. Coli. |
| <b>Preservatives</b>       | Purified antibody in PBS containing 0.03% sodium azide.          |
| <b>Form/Appearance</b>     | Purified antibody in PBS containing 0.03% sodium azide.          |
| <b>Concentration</b>       | 1mg/ml                                                           |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.      |
| <b>Note</b>                | For research use only                                            |
| <b>Isotype</b>             | Mouse IgG1                                                       |
| <b>Clonality</b>           | Monoclonal                                                       |
| <b>Entrez</b>              | <a href="#">4151</a>                                             |
| <b>Dilution Range</b>      | ELISA: 1/10000, IHC-P: 1/200 - 1/1000, WB: 1/500 - 1/2000        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |



Western blot analysis of truncated Myogl...



Immunohistochemical analysis of paraffin...