



Product Datasheet

GFP Antibody (orb69018)


Description	Mouse monoclonal antibody to GFP
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Immunogen	Purified recombinant fragment of GFP expressed in E. Coli.
Preservatives	Purified antibody in PBS containing 0.03% sodium azide.
Form/Appearance	Purified antibody in PBS containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	This product will ship in a box containing blue ice at a temperature of 4°C.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	27kDa
Dilution Range	ELISA: 1/10000, IHC-P: 1/200 - 1/1000, FACS: 1/200 - 1/400, WB: 1/500 - 1/2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from HCC827 cells...



Flow cytometric analysis of HCC827 cells...



Immunocytochemistry analysis of HCC827 cells...