



Product Datasheet

LCN1 Antibody (orb69147)


Description	Mouse monoclonal antibody to LCN1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC
Immunogen	Purified recombinant fragment of LCN1 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	19.2kDa
Entrez	3933
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Western blot analysis of HEK293 (lane1) ...



Immunohistochemical analysis of paraffin...