



## Product Datasheet

### ISL1 Antibody (orb69159)

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to ISL1                                 |
| <b>Species/Host</b>        | Mouse                                                             |
| <b>Reactivity</b>          | Human                                                             |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Tested Applications</b> | ELISA, ICC                                                        |
| <b>Immunogen</b>           | Purified recombinant fragment of human ISL1 expressed in E. Coli. |
| <b>Preservatives</b>       | Ascitic fluid containing 0.03% sodium azide.                      |
| <b>Concentration</b>       | 1mg/ml                                                            |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.       |
| <b>Note</b>                | For research use only                                             |
| <b>Isotype</b>             | Mouse IgG1                                                        |
| <b>Clonality</b>           | Monoclonal                                                        |
| <b>MW</b>                  | 39kDa                                                             |
| <b>Entrez</b>              | <a href="#">3670</a>                                              |
| <b>Dilution Range</b>      | ICC: 1/200 - 1/1000; ELISA: 1/10000                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |



Confocal immunofluorescence analysis of ...



Western blot analysis of full-length ISL...