



Product Datasheet

C5aR2 antibody (orb691931)

Description	Mouse monoclonal antibody to C5aR2
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	L1.2 cells transfected with human C5aR2
Target	C5aR2
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-6 µg/ml
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9P296
Entrez	27202
Dilution Range	Flow cytometry: Recommended dilution: 2-6 µg/ml
Expiration Date	12 months from date of receipt.



Separation of monocytes stained using an...



Flow cytometry surface staining pattern ...