

Product Datasheet

Novel Coronavirus Nucleoprotein (SARS-CoV-2 N) ELISA Kit (orb692157)

Description	Novel Coronavirus Nucleoprotein (SARS-CoV-2 N) ELISA Kit
Range	3.12 ng/mL-200 ng/mL
Target	SARS-CoV-2 N
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, swabs
Uniprot ID	P0DTC9
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.