

Product Datasheet

METTL3 Antibody (orb692196)

Description
METTL3 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IF, IHC, WB

Applications
Immunogen

E.coli-derived human METTL3 recombinant protein (Position: E71-L580).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


64 kDa

Uniprot ID
Q86U44
Dilution Range


Western blot, 0.25-0.5µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Human, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 5µg/ml, Human
Flow Cytometry, 1-3µg/1x10⁶ cells, Human Direct ELISA, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.



Western blot analysis of METTL3 using anti-METTL3...



IHC analysis of METTL3 using anti-METTL3...




Figure 3 IHC analysis of METTL3 using anti-METTL3...