

Product Datasheet

PROM1/CD133 Antibody (monoclonal, 7F13) (orb692236)

Description

PROM1/CD133 Antibody (monoclonal, 7F13)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

E.coli-derived human PROM1/CD133 recombinant protein (Position: P531-H865).
HumanPROM1/CD133 shares 61% amino acid (aa) sequence identity with mouse PROM1/CD133.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt.
After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

120 kDa

Uniprot ID


[O43490](#)

Dilution Range

Western blot, 0.25-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.



Western blot analysis of PROM1/CD133 protein.