

Product Datasheet

HighFidelity GREEN PCR Labeling Testkit (orb692240)

Description	HighFidelity GREEN PCR Labeling Testkit is designed to produce randomly Fluorescein-, ATTO488-and AF...
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: Fluorescein: λ_{exc} 492 nm, λ_{em} 517 nm, ϵ 83.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5), ATTO488: λ_{exc} 500 nm, λ_{em} 520 nm, ϵ 90.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5), AF488: λ_{exc} 494 nm, λ_{em} 515 nm, ϵ 73.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.