

Product Datasheet

HighFidelity YELLOW PCR Labeling Testkit (orb692244)

Description	HighFidelity YELLOW PCR Labeling Testkit is designed to produce randomly Cy3-and AF555-modified DNA ...
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: Cy3: λ_{exc} 550 nm, λ_{em} 570 nm, ϵ 150.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5), AF555: λ_{exc} 555 nm, λ_{em} 572 nm, ϵ 155.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.