


Product Datasheet

PRDM1 Antibody (orb69241)

Description	Mouse monoclonal antibody to PRDM1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant fragment of human PRDM1 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	88kDa
Entrez	639
Dilution Range	WB: 1/500 - 1/2000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of Raji (lane1)