



Product Datasheet

MUM1 Antibody (orb69264)


Description	Mouse monoclonal antibody to MUM1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Purified recombinant fragment of human MUM1 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	78.6kDa
Entrez	84939
Dilution Range	IHC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of human MUM1 reco...



Immunohistochemical analysis of paraffin...