


## Product Datasheet

### CCND1 Antibody (orb69273)

|                            |                                                                    |
|----------------------------|--------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CCND1                                 |
| <b>Species/Host</b>        | Mouse                                                              |
| <b>Reactivity</b>          | Human                                                              |
| <b>Conjugation</b>         | Unconjugated                                                       |
| <b>Tested Applications</b> | ELISA                                                              |
| <b>Immunogen</b>           | Purified recombinant fragment of human CCND1 expressed in E. Coli. |
| <b>Preservatives</b>       | Ascitic fluid containing 0.03% sodium azide.                       |
| <b>Concentration</b>       | 1mg/ml                                                             |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.        |
| <b>Note</b>                | For research use only                                              |
| <b>Isotype</b>             | Mouse IgG1                                                         |
| <b>Clonality</b>           | Monoclonal                                                         |
| <b>MW</b>                  | 33.7kDa                                                            |
| <b>Entrez</b>              | <a href="#">595</a>                                                |
| <b>Dilution Range</b>      | ELISA: 1/10000                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                    |



Western blot analysis of CCND1-hlgFc tr...