

Product Datasheet

Bovine Plg ELISA Kit (orb692733)

Description	Sandwich Bovine Plasminogen (Plg) ELISA Kit
Reactivity	Bovine
Range	7.81-500ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P06868
Sensitivity	3.2ng/mL
Entrez	280897
Expiration Date	Please enquire.