

Product Datasheet

Human LAMA1 ELISA Kit (orb692752)

Description	Sandwich Human Laminin Subunit Alpha-1 (LAMA1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P25391
Sensitivity	31pg/mL
Entrez	284217
Expiration Date	Please enquire.