

Product Datasheet

Human MC4R, ELISA, Kit (orb692756)

Description	Sandwich Human Melanocortin 4 Receptor (MC4R) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P32245
Sensitivity	0.060ng/mL
Expiration Date	Please enquire.