

Product Datasheet

Human MC4R, ELISA, Kit (Ready to Use) (orb692757)

Description	Sandwich Human Melanocortin 4 Receptor (MC4R) ELISA Kit Ready to Use: Pre-diluted detecti...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P32245
Sensitivity	0.060ng/mL
Expiration Date	Please enquire.