

Product Datasheet

Mouse CHGA ELISA Kit (orb692765)

Description	Sandwich Mouse Chromogranin A (CHGA) ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P26339
Sensitivity	22.3pg/mL
Entrez	12652
Expiration Date	Please enquire.