

Product Datasheet

Mouse CHGA ELISA Kit (Ready to Use) (orb692766)

Description	Sandwich Mouse Chromogranin A (CHGA) ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P26339
Sensitivity	22.3pg/mL
Entrez	12652
Expiration Date	Please enquire.