

Product Datasheet

Mouse CD8a ELISA Kit (orb692767)

Description	Sandwich Mouse Cluster Of Differentiation 8a (CD8a) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P01731
Sensitivity	32.6pg/mL
Entrez	12525
Expiration Date	Please enquire.