

Product Datasheet

Mouse EDN2 ELISA Kit (orb692771)

Description	Sandwich Mouse Endothelin 2 (EDN2) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P22389
Sensitivity	11.2pg/mL
Expiration Date	Please enquire.