

Product Datasheet

Mouse HSP90aB1 ELISA Kit (orb692777)

Description	Sandwich Mouse Heat Shock Protein 90kDa Alpha B1 (HSP90aB1) ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P11499
Sensitivity	0.56ng/mL
Expiration Date	Please enquire.