

Product Datasheet

Porcine F5 ELISA Kit (orb692792)

Description	Sandwich Porcine Coagulation Factor V (F5) ELISA Kit
Reactivity	Porcine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	plasma and other biological fluids.
Uniprot ID	Q9GLP1
Sensitivity	0.120ng/mL
Entrez	397217
Expiration Date	Please enquire.