

Product Datasheet

Rat GAP43 ELISA Kit (orb692800)

Description	Sandwich Rat Growth Associated Protein 43 (GAP43) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P07936
Sensitivity	0.057ng/mL
Entrez	29423
Expiration Date	Please enquire.