

Product Datasheet

Rat GAP43 ELISA Kit (Ready to Use) (orb692801)

Description	Sandwich Rat Growth Associated Protein 43 (GAP43) ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P07936
Sensitivity	0.057ng/mL
Entrez	29423
Expiration Date	Please enquire.