

Product Datasheet

Bovine IL12A ELISA Kit (orb692808)

Description	Sandwich Bovine Interleukin 12A (IL12A) ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	5.9pg/mL
Expiration Date	Please enquire.