

Product Datasheet

Bovine IL12A ELISA Kit (Ready to Use) (orb692809)

Description	Sandwich Bovine Interleukin 12A (IL12A) ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	5.9pg/mL
Expiration Date	Please enquire.