

## Product Datasheet

### Goat IL2 ELISA Kit (orb692838)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>       | Sandwich Goat Interleukin 2 (IL2) ELISA Kit                                              |
| <b>Reactivity</b>        | Goat                                                                                     |
| <b>Range</b>             | 15.6-1000pg/mL                                                                           |
| <b>Form/Appearance</b>   | Traditional                                                                              |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                    |
| <b>Application notes</b> | serum, plasma, tissue homogenates and other biological fluids.                           |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates and other biological fluids.                           |
| <b>Sensitivity</b>       | 6.9pg/mL                                                                                 |
| <b>Expiration Date</b>   | Please enquire.                                                                          |