

Product Datasheet

Human aEP ELISA Kit (orb692845)

Description	Competitive Inhibition Human Alpha-Endorphin (aEP) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluid.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Sensitivity	4.45pg/mL
Expiration Date	Please enquire.