

Product Datasheet

Human CIRBP ELISA Kit (orb692852)

Description	Sandwich Human Cold-Inducible RNA-Binding Protein (CIRBP) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q14011
Sensitivity	32.6pg/mL
Entrez	1153
Expiration Date	Please enquire.