

Product Datasheet

Human dSIP ELISA Kit (orb692859)

Description	Sandwich Human Delta-Sleep Inducing Peptide (dSIP) ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	47pg/mL
Expiration Date	Please enquire.