

Product Datasheet

Human Human Macrosialin ELISA Kit (orb692869)

Description	Sandwich Human Macrosialin ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P34810
Sensitivity	0.120ng/mL
Entrez	968
Expiration Date	Please enquire.