

Product Datasheet

Human PANX1 ELISA Kit (orb692876)

Description	Sandwich Human Pannexin-1 (PANX1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q96RD7
Sensitivity	0.060ng/mL
Entrez	24145
Expiration Date	Please enquire.