

Product Datasheet

Human PRSS3 ELISA Kit (Ready to Use) (orb692883)

Description	Sandwich Human Protease, Serine 3 (PRSS3) ELISA KitReady to Use: Pre-diluted detection reagents for ...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	26.9pg/mL
Expiration Date	Please enquire.