

Product Datasheet

Mouse CASP12 ELISA Kit (orb692896)

Description	Sandwich Mouse Caspase 12 (CASP12) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	008736
Sensitivity	0.065ng/mL
Entrez	12364
Expiration Date	Please enquire.