

Product Datasheet

Mouse C1qA ELISA Kit (orb692898)

Description	Sandwich Mouse Complement 1q Subcomponent Subunit A (C1qA) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P98086
Sensitivity	0.123ng/mL
Entrez	12259
Expiration Date	Please enquire.