

Product Datasheet

Mouse GRIA1 ELISA Kit (Ready to Use) (orb692907)

Description	Sandwich Mouse Glutamate Receptor, Ionotropic, AMPA 1 (GRIA1) ELISA KitReady to Use: Pre-diluted det...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P23818
Sensitivity	0.051ng/mL
Entrez	14799
Expiration Date	Please enquire.