

Product Datasheet

Mouse GP4 ELISA Kit (orb692917)

Description	Sandwich Mouse Platelet Membrane Glycoprotein IV (GP4) ELISA Kit
Reactivity	Mouse
Range	0.39-25ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q08857
Sensitivity	0.18ng/mL
Entrez	12491
Expiration Date	Please enquire.