

Product Datasheet

Mouse TNFSF12 ELISA Kit (orb692922)

Description	Sandwich Mouse Tumor Necrosis Factor Ligand Superfamily, Member 12 (TNFSF12) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O54907
Sensitivity	31.6pg/mL
Entrez	21944
Expiration Date	Please enquire.