

Product Datasheet

Rat C4BPα ELISA Kit (orb692927)

Description	Sandwich Rat C4 Binding Protein Alpha (C4BPα) ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q63514
Sensitivity	25.6pg/mL
Entrez	24235
Expiration Date	Please enquire.