

Product Datasheet

Rat LIPF ELISA Kit (orb692931)

Description	Sandwich Rat Lipase, Gastric (LIPF) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P04634
Sensitivity	32.6pg/mL
Entrez	50682
Expiration Date	Please enquire.