

Product Datasheet

Rat S100 ELISA Kit (orb692933)

Description	Sandwich Rat S100A1 Calcium Binding Protein (S100) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P35467
Sensitivity	0.57ng/mL
Entrez	295214
Expiration Date	Please enquire.